<b>DE</b>	<b>EMSER® INHALATOR PRO</b>	
	Druckluftvernebler .....	2
<b>FR</b>	<b>EMSER® INHALATEUR PRO</b>	
	Nébuliseur à air comprimé .....	8
<b>IT</b>	<b>EMSER® INALATORE PRO</b>	
	Nebulizzatore ad aria compressa..	14
<b>EN</b>	<b>EMSER® INHALATOR PRO</b>	
	Compressed air atomiser .....	20

Abbildung 1: Kompressor



Abbildung 2: Zubehör

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Globalcare Medical Technology Co., Ltd.,  
7<sup>th</sup> Building, 39 Middle Industrial Main Road,  
European Industrial Zone, Xiaolan Town,  
528415 Zhongshan City, Guangdong Province, PRC  
Tel.: +86 760 22589901  
<http://www.globalcare.com.hk/contact/>

Donawa Lifescience Consulting Srl, Piazza  
Albania, 10, 00153 Rome, Italia  
<https://www.donawa.com/wli/main/contatti.index>

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

**Figure 1: Compressor**

**Figure 2: Accessories**

1. Compressor
  - 1a Operating switch
  - 1b Compressed air hose connection point
  - 1c Air filter compartment
2. Compressed air hose with straight or angled end
3. Atomiser set
- 3a: Nebuliser head
4. Mouthpiece
5. Adult face mask
6. Children's face mask
7. Nosepiece
8. Replacement air filter
9. Connecting the compressed air hose with the atomiser set
10. Assembled atomiser set
11. Connecting the face mask to the atomiser set

## Dear customer,

The EMSER® Inhalator Pro is an aerosol therapy system for domestic use. This device is suited for atomising liquids and liquid medication (aerosols) and for the treatment of the upper and lower airways. By atomising and inhaling the inhalation solution recommended by your doctor, you can prevent diseases affecting the airways, or in the case that you contract such an illness, you can alleviate symptoms and speed up your recovery. Ask your doctor or pharmacist for further information about the potential applications.


Read these instructions carefully to understand all functions and safety instructions. If you have any questions or problems or wish to order spare parts, please contact our customer service (see Section 7).

### Intended purpose:

The Emser® Inhalator Pro is an aerosol therapy system for the treatment of respiratory diseases for domestic use.

The Emser® Inhalator Pro is used to produce medical aerosols using a compressed air operated nebulizer kit.

## 1. Explanation of symbols

	Medical device
	Electronic devices must be disposed of in accordance with the applicable regulations and may not be put in the bin
	Observe the instructions for use
	Important
	Type BF application part
	Protection class II device
	Serial number
	Item number
	Manufacturer
	Distributor
	EC Authorised Representative
	ON

	OFF
<b>IP21</b>	Degree of protection against ingress of solid foreign bodies and damage caused by infiltration of water
<b>CE 0123</b>	This product satisfies the requirements of the applicable European and national directives.
	Storage and transport conditions
	Operating conditions
	Dispose of packaging in an environmentally friendly manner
	Single patient – multiple use
	Importer
	Authorized representative in Switzerland

## 2. Important safety instructions

- The EMSER® Inhalator Pro is only intended for use as described in these instructions. The manufacturer is not liable for damage caused by improper use.
- Take good care of these instructions for future reference.
- The EMSER® Inhalator Pro is not suitable for use in the presence of flammable anaesthetic mixtures with air, oxygen or nitrogen oxide.
- The EMSER® Inhalator Pro should only be used with the original accessories shown in these instructions.
- Do not use the EMSER® Inhalator Pro if you suspect damage or notice anything unusual.
- The EMSER® Inhalator Pro must not be opened.
- The EMSER® Inhalator Pro contains sensitive components and must be handled with care. See the storage and operating conditions in the section “Technical data”.
- Protect the EMSER® Inhalator Pro from:
  - water and moisture
  - extreme temperatures
  - knocks and falling
  - dirt and dust
  - direct sunlight
- Adhere to the safety regulations for electrical devices, in particular:
  - Do not touch the EMSER® Inhalator Pro with wet or damp hands.
  - Place the EMSER® Inhalator Pro on a stable and horizontal surface for operation.
  - Do not pull on the power cable or the EMSER® Inhalator Pro itself to disconnect the mains plug.

- The EMSER® Inhalator Pro is disconnected from the power supply by unplugging the mains plug. Ensure that it is accessible during operation.
- Before inserting the EMSER® Inhalator Pro, check whether the nominal voltage on the type plate on the underside of the EMSER® Inhalator Pro corresponds to the mains voltage.
- If the supplied mains plug does not fit into the socket, ask a qualified person to replace the plug with a suitable one. As a general rule, the use of single or multi plug adapters or extension cables is not recommended. If their use is indispensable, models that correspond to the safety regulations must be used. The maximum voltage stated on the adapters and extension cables must not be exceeded.
- Do not leave the EMSER® Inhalator Pro plugged in when not in use. Disconnect it from the mains.
- The installation must be carried out in accordance with the manufacturer's specifications. Improper installation can result in harm to people, animals or objects, for which the manufacturer is not liable.
- Do not replace the EMSER® Inhalator Pro power cable. If the power cable is damaged, contact a specialist company authorised by the manufacturer to provide a replacement.
- To avoid overheating, the power cable must always be completely unwound.
- Switch off the device and unplug from the mains before servicing or cleaning it.
- The device is maintenance-free.
- Only use the medications prescribed by your doctor and follow their instructions regarding dosage, duration and frequency of treatment.
- Only use the treatment recommended by your doctor for your illness.
- Do not use the nosepiece unless specifically instructed to do so by your doctor and DO NOT insert the nozzles into your nose. Instead, place them as close as possible to your nose.
- Position the EMSER® Inhalator Pro so that it is easily accessible and you can easily pull out the plug.
- The atomiser set (3) and other accessories are for single patient reuse. The device is for multi-patient reuse, please follow the instructions in section 4 "Cleaning and disinfection".
- The atomiser must not be tilted above 60 degrees during use.
- Do not use the EMSER® Inhalator Pro in the vicinity of strong electromagnetic fields, e.g. those generated by mobile phones and radio devices. When operating the atomiser, keep a minimum distance of 0.3 m from those devices.

#### Contraindications:

- No contraindications exist to the administration of aerosols by inhalation. Contraindications related to the medicine used must be checked with the medicines package leaflet.

- The EMSER® Inhalator Pro is not indicated to use with quick-relief medications during life-threatening asthma attack.
- The EMSER® Inhalator Pro is designed to nebulize solution and suspension liquids. Nebulizing system is not suitable for use in anaesthetic breathing system or a ventilator breathing system.

Consult your physician in case of doubts.

Do not allow children to use the EMSER® Inhalator Pro unattended; some parts are small enough to be swallowed. Be aware of the risk of strangulation, as the device has cables and tubes.

Using this device is not a substitute for medical consultation and treatment. Always consult your doctor immediately in the event of pain or illness. If you have health concerns of any kind, consult your GP.

The EMSER® Inhalator Pro is not a toy! Keep it out of reach of children and do not leave children alone during treatment.

Any serious incident that has occurred in relation to the product should be reported to the manufacturer and the competent authority of the Member State in which the user is established or to the European Authorised Representative (EC REP). National vigilance contact points: <https://ec.europa.eu/growth/sectors/medical-devices/contacts/>

### 3. Preparation and use of the EMSER® Inhalator Pro

Before using the EMSER® Inhalator Pro for the first time, we recommend cleaning it as described in the "Cleaning and disinfection" section.

1. Unpack the device and check that all parts are present.
2. If not already assembled on delivery, assemble the atomiser set (3) and connect the compressed air hose (2) with the straight end to the atomiser set (3, 10) (see Figure 9)
3. Carefully connect the angled end of the compressed air hose (2) to the connection point on the compressor (1b).
4. Fill the atomiser (3) with an inhalation solution recommended by your doctor or pharmacist. Do not exceed the maximum fill level of 8 ml indicated by the "max" marking.
5. Connect the atomiser set to the desired accessory, mouthpiece (4), face mask (5, 6, see Figure 11) or nosepiece (7).
  - The mouthpiece provides better delivery of the medication to the lungs.

- Choose either the adult face mask (5) or children's face mask (6) and make sure that the mouth and nose are completely covered.
  - Use all accessories, including nosepiece (7), as prescribed by your physician.
6. Insert the power cable into the socket.
  7. Turn the operating switch (1a) on to start the treatment.
  8. During inhalation, sit upright and relaxed at a table, not in an armchair. This is to avoid compressing the airways and compromising the effectiveness of the treatment. **Do not inhale while lying down.** Stop inhalation if you feel uncomfortable.
  9. After the period of inhalation recommended by your physician, turn the EMSER® Inhalator Pro off with the operating switch (1a) and remove the mains plug.
  10. Empty the atomiser (3) of any remaining medication and clean it as described in the section "Cleaning and disinfection".

 **The EMSER® Inhalator Pro is designed for interval operation with 30 minutes operation / 30 minutes pause.** Switch off the device after 30 minutes at the latest, and wait 30 minutes before continuing treatment.

 This device does not need to be calibrated.

 The device may not be modified.

 Essential oils, cough syrups, gargling solutions and drops to be used as a rub or in a steam bath are wholly unsuitable for inhalation using a nebuliser. These additives are often viscous and can impair the correct functioning of the device and therefore the effectiveness of the application in the long term. For individuals with a hypersensitive bronchial system, medications containing essential oils may under certain conditions cause an acute bronchospasm (a sudden cramp-like restriction of the bronchi with shortness of breath). Consult your doctor or pharmacist in relation to this matter.

## 4. Cleaning and disinfection

After treatment and if not used for the longer term, thoroughly clean all components as described below to remove any drug residue and possible contamination. Use a soft dry cloth and, if necessary, a non-abrasive cleaning product to clean the compressor.

 Ensure that the inner parts of the EMSER® Inhalator Pro do not come into contact with liquid, and remove the plug from the mains before cleaning.

### Cleaning and disinfection of accessories

Carefully follow the cleaning and disinfection instructions for the accessories, as they are essential for the performance of the device and the success of the treatment.

### Before first use or after each further treatment,

disassemble the atomiser (3) into its individual parts by turning the upper part to the left and detaching it from the lower part, then remove the cover and the cone. Wash the disassembled atomiser components, the mouthpiece (4) and the nosepiece (7) with tap water and immerse them in boiling water for 5 minutes, avoiding contact between the accessories and the hot pan base. Reassemble the atomiser components, connect to the compressed air hose with the device switched on and let it run for 10 to 15 minutes.

For disinfection of the masks (5, 6), wash them with warm water and immerse the masks in 70 % ethanol for 10 minutes. Afterwards, wash the masks with tap water again.

### Storage

#### Drying

Place the individual parts on a clean, dry and absorbent surface and leave them to dry completely (at least 4 hours).

#### Note

Please be aware that the parts must be completely dried after cleaning, as otherwise the risk of bacterial growth is increased.

For storing the device, put the parts together again once these are completely dry. Place the parts in a dry, sealed container.

 Wash the masks with warm water.

 The compressed air hose and masks must not be boiled or autoclaved (steam sterilisation process between 110° C and 140° C).

## 5. Maintenance, care and service

**Replacement parts and wearing parts can be purchased from the service address indicated in Section 7**

Designation	REF
Mouthpiece (PP, silicone), adult mask (PVC, aluminium), children's mask (PVC, aluminium), atomiser (PP, silicone), compressed air hose (PVC), nosepiece (PP, silicone), filters (polyester cotton)	533558

### Replacing the atomiser

Replace the atomiser (3) if it is deformed or damaged or if the nebuliser head (3a) is blocked by dried medication, dust, etc. following a long break from treatment. The atomiser should be replaced after a period of half a year to one year, depending on the frequency of treatment.

 Only use genuine atomisers!

## Replacing the air filter

In normal operating conditions, the air filter must be replaced after approx. 200 operating hours or one year. Check the air filter regularly (after about 10 to 12 treatments). If it is grey or brown or wet, remove the filter and replace it. To do this, carefully loosen the cover of the air filter compartment (1c) and replace the filter (8). Close the lid (1c) again.

☞ The filter must not be cleaned and reused.

☞ The filter must not be removed for servicing while the EMSER® Inhalator Pro is in use.

 Only use genuine filters! Do not use the device without a filter!

## 6. Malfunctions and countermeasures

### The device will not switch on.

- Ensure that the power cable is inserted correctly.
- Ensure that the operating switch (1a) is in position I.
- Ensure that the unit has been used within the limits of the operating time specified in these instructions.

### The atomiser is working poorly or not working at all

- Ensure that the compressed air hose (2) is correctly connected at both ends.
- Check that the compressed air hose is not crushed, kinked, dirty or blocked. Replace it with a new one if necessary.
- Check that the atomiser (3) is fully assembled and that the coloured nebuliser head (3a) is attached correctly and not blocked.
- Ensure that the medication has been added to the device.

## 7. Warranty/service

We offer a 5-year product warranty from the date of purchase against material and production faults. The warranty does not cover:

- Damage due to improper operation.
- Wearing parts.
- Cases of personal negligence on the part of the customer.
- Cases in which the device has been opened by a non-authorised workshop.

The warranty does not affect the customer's statutory rights. To assert a warranty claim within the warranty period, the customer must provide proof of purchase. The warranty claim must be submitted within a period of 5 years from the date of purchase to Sidroga GfGmbH, Bad Ems, for Germany or Austria, or to Sidroga AG, Rheinfelden, for Switzerland.

In case of complaints, please contact our service department using the following contact details:

## Germany/Austria:

Service hotline:

Tel.: +49 (0)2603/9604-710

Fax: +49 (0)2603/9604-711

E-Mail: info@emser.de

www.emser.de

If requested, please send the faulty product to the following address:

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
E-Mail: apothekenservice@sidroga-pharma.com

www.emser.ch

If requested, please send the faulty product to the following address:

Sidroga AG  
Weidenweg 15  
CH-4310 Rheinfelden  
SWITZERLAND

## 8. Technical data

Type: GCE843

**AEROSOL PERFORMANCES**

**ACCORDING TO**

**EN27427:2019 based on**

**adult ventilatory pattern**

**with Sodium Fluoride (NaF):**

- 1) Aerosol output: 0.11 ml
- 2) Aerosol output rate: 0.07 ml/min
- 3) Percentage of fill volume emitted per min: 3.5 %
- 4) Residual Volume: 1.68 ml

- 5) Particle size (MMAD): 3.07 µm
- 6) GSD (geometric standard deviation): 2.03
- 7) RF (respirable fraction < 5 µm): 75 %
- 8) Large particle range (> 5 µm): 25 %
- 9) Middle particle range (2 to 5 µm): 47.8 %
- 10) Small particle range (< 2 µm): 27.2 %

Donawa Lifescience Consulting Srl, Piazza Albania, 10, 00153 Rome, Italy  
<https://www.donawa.com/wli/main/contatti.index>

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Sidroga AG, Weidenweg 15, CH-4310 Rheinfelden

Date of information / latest revision: 2022-02

**Atomisation rate:** 0.35 ml/min  
**Maximum pressure :** 2.2 bar  
**Noise emission:** 52 dB(A)  
**Power source:** AC 230 V 50 Hz, ≤ 1000 mA  
**Atomisation capacity:** min. 2 ml; max. 8 ml  
**Operating time limits:** 30 minutes on/ 30 minutes off  
**Operating conditions:** 10–40 °C  
 10–95 % relative maximum humidity  
 700–1060 hPa

**Storage and shipping conditions:** -25 – +70 °C  
 10–95 % relative maximum humidity  
 700–1060 hPa

**Weight:** 1500 g  
**Dimensions:** 150×150×95 mm  
**IP class:** IP21  
**Relevant standards:** EN 13544-1; EN, 60601-1; EN60601-2; EN 60601-1-6; IEC 60601-1-11

**Expected operating hours:** 1,000 hours

This device complies with Medical Devices Directive 93/42/EEC.

Protection against electric shock: Class II atomiser, mouthpiece and masks are type BF application parts.

Technical data is subject to change without notice.

## 9. Disposal

### General disposal

For environmental reasons, do not dispose of the device in your household waste at the end of its service life. Dispose of the device at appropriate collection points in your country. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Globalcare Medical Technology Co., Ltd.,  
 7<sup>th</sup> Building, 39 Middle Industrial Main Road,  
 European Industrial Zone, Xiaolan Town,  
 528415 Zhongshan City, Guangdong Province, PRC  
 phone: +86 760 22589901  
<http://www.globalcare.com.hk/contact/>

## AC Compressor Atomiser – Guidance and manufacturer’s declaration

**Table 1 – Electromagnetic emissions**

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

## AC Compressor Atomiser – Guidance and manufacturer’s declaration

**Table 2 – Electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below.

The customer or the user of the device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	Power supply lines ±2 kV input/output lines ±1 kV	Power supply lines ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s) ±1 kV. line(s) to earth ±2 kV. 100 kHz repetition frequency	line(s) to line(s) ±1 kV. 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycle  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycle  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5s	5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycle  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycle  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5s	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: $U_T$ is the a.c. mains voltage prior to application of the test level.			

## AC Compressor Atomiser – Guidance and manufacturer’s declaration

**Table 3 – Electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below.

The customer or the user of the device should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance	
Conducted RF IEC61000-4-6	150 KHz to 80 MHz: 3 Vrms 6 Vrms (in ISM and amateur radio bands) 80 % Am at 1 kHz	150 KHz to 80 MHz: 3 Vrms 6 Vrms (in ISM and amateur radio bands) 80 % Am at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device ,including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. <b>Recommended separation distances:</b> $d=1.2\sqrt{P}$ ; $d=2\sqrt{P}$	
Radiated RF IEC61000-4-3	10 V/m	10 V/m	80 MHz to 800 MHz: $d=1.2\sqrt{P}$ 800 MHz to 2.7 GHz. $d=2.3\sqrt{P}$	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## AC Compressor Atomiser – Guidance and manufacturer’s declaration

**Table 4 – Recommended separation distances between portable and mobile RF communications equipment and the device**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d=1.2\sqrt{P}$	150 kHz to 80 MHz (in ISM and amateur radio bands) $d=2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.7GHz $d=2.3\sqrt{P}$
0.01	0.12	0.2	0.12	0.23
0.1	0.38	0.632	0.38	0.73
1	1.2	2	1.2	2.3
10	3.8	6.32	3.8	7.3
100	12	20	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## AC Compressor Atomiser – Guidance and manufacturer’s declaration

**Table 5 – Electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below.

The customer or the user of the device should assure that it is used in such an environment.

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380 – 390	TETRA 400	Pulse modulation b) 18 Hz	1.8	0.3	27
	450	380 – 390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0.2	0.3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0.3	28
	870						
	930						
	1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28
	5240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
	5240						
	5785						

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = \frac{6}{d} \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.